ClinCal[®] - Calibrator

Urine Calibrator lyophilised / Urin-Kalibrator lyophilisiert

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Intended use:

ClinCal[®] Urine Calibrators are used for calibration of your analyzing system in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. This lyophilised calibrator is based on human urine and is available in a range of concentration that allows a comfortable one-point-calibration. After reconstitution the calibrator has to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 5.0 ml of analytical grade or double distilled water to the vial and close it well. Then dissolve the lyophilised material completely. When all material is dissolved, the solution is ready to use.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C.

After reconstitution the stability of the analytes is:

- 8 hours when stored at 15 - 30 °C
- 7 days when stored at 2 - 8 °C
- 30 days when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinCal[®] Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.61 %.

Mean values:

The mean values have been established in additional independent reference laboratories with supervision of RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. For the determination of all parameters, chromatographic techniques were performed. Even if you use different techniques of analysis, you should find your control results within the ranges indicated on the data sheet, after having calibrated your analyzing system with this ClinCal[®] urine calibrator.

Pack size:

ClinCal[®] Urine Calibrator
5 x 5 ml, **order no.: 9922**

Precautions:

This calibrator was manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH