ClinCal® - Calibrator

Whole Blood Calibrator Lyophilised / Vollblut-Kalibrator lyophilisiert

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Intended use:

ClinCal® Whole Blood Calibrators are used for calibration of your analyzing system in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. This lyophilised calibrator is based on human whole blood and is available in a range of concentration, which allows a comfortable one-point-calibration. After reconstitution the calibrator has to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 5.0 ml of analytical grade or double distilled water to the vial and close it well. Then dissolve the lyophilised material completely. For this purpose we recommend a so-called roller mixer; a treatment of 30 min with such an instrument is sufficient. Otherwise incubate the samples for 1 to 2 hours at room temperature by shaking it carefully from time to time. Prior to use mix the calibrator well again.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After reconstitution the stability of the analytes is:

- 8 hours when stored at 15 - 30 °C
- 4 days when stored at 2 - 8 °C
- 30 days when stored below -18 °C (avoid repeated freezing and thawing)

If performing mercury analyses, use only freshly dissolved material.

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.39 %.

Mean values:

The mean values have been established in independent reference laboratories with supervision of RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. For the determination of the trace elements, atomic absorption spectrometry (AAS) was used as well as ICP-MS. Even if you use different techniques of analysis, you should find your control results within the ranges indicated on the data sheet, after having calibrated your analyzing system with this ClinCal® Whole Blood Calibrator.

Pack size:

ClinCal® Whole Blood Calibrator
5 x 5 ml, **order no.: 9943**

Precautions:

The human whole blood which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV 1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV1- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the calibrator should be considered as potentially infectious and treated with appropriate care.

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*: For Iodide we could not establish any mean values, due to the low number and the high imprecision of the available results. /

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In Germany Calcium, Magnesium, Potassium and Sodium are subject to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK). For assessment of the internal quality control, please also refer to table B 1a of this guideline.

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CHEMICALS + INSTRUMENTS GmbH