ClinCal® - Calibrator Serum Calibrator Set lyophilised / Serum-Kalibrator-Set lyophilisiert

FOR BENZODIAZEPINES / FÜR BENZODIAZEPINE

Intended use:

ClinCal® Serum Calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on human serum and are available as a calibrator set with 4 levels, including a blank. After reconstitution the calibrators have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 1.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C.

After reconstitution the stability of the analytes is:

- 1 day when stored at 15 - 30 °C (1 hour for Zopiclone)
- 7 days when stored at 2 - 8 °C (3 hours for Zopiclone)
- 90 days when stored below -18 °C (also for Zopiclone) (avoid repeated freezing and thawing)

Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.83 %.

Mean values:

The mean values have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinCal® Serum Calibrator Set
Level 0 - 3
4 x 1 x 1 ml, **order no.: MS6013**

Precautions:

The human serum which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV- and HCV-RNA, HBV-DNA. Nevertheless, the serum calibrators should be considered as potentially infectious and treated with appropriate care.

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2019-10

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Level 0

This blank-level is manufactured from the same pool of human serum like the others, but it was not spiked with any Benzodiazepine. Therefore the concentration of each of the 35 analytes is zero. /

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Method of Analysis / Analysenmethode: LC-MS/MS

RECIPE CHEMICALS + INSTRUMENTS GmbH
Dessauerstr. 3 D-80992 Munich / Germany
Tel.: +49 / 89 / 54 70 81 - 0 Fax: +49 / 89 / 54 70 81 - 11
Internet: www.recipe.de e-mail: info@recipe.de