ClinCal® - Calibrator
Serum Calibrator Set lyophilised / Serum-Kalibrator-Set lyophilisiert

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Intended use:

ClinCal® Serum Calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on human serum and are available as a calibrator set with 4 levels, including a blank. After reconstitution the calibrators have to be prepared like patient samples in one series of analyses.

After reconstitution the stability of the analytes is:

- 1 day when stored at 15 - 30 °C
- 7 days when stored at 2 - 8 °C
- 90 days when stored below -18 °C (avoid repeated freezing and thawing)

Mean values:

The mean values have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinCal® Serum Calibrator Set
 Level 0 - 3
 4 x 1 x 1 ml, **order no.: MS8113**

Reconstitution:

Add exactly 1.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 1.15 %.

Precautions:

The human serum which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen. Nevertheless, the serum calibrators should be considered as potentially infectious and treated with appropriate care.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C.

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2021-05

Level 0

This blank-level was manufactured from the same pool of serum like Level 1 - 3. It was not spiked with any analytes. Level 0 shall not be included in the calibration curve, rather shall be set as an „unknown“.

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Technique of analysis / Analysentechnik: LC-MS/MS

RECIPE CHEMICALS + INSTRUMENTS GmbH
Dessauerstr. 3 D-80992 Munich / Germany
Tel.: +49 / 89 / 54 70 81 - 0 Fax: +49 / 89 / 54 70 81 - 11
Internet: www.recipe.de e-mail: info@recipe.de