ClinChek® - Control Hemolysate Control / Kontrollhämolysat

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Intended use:

ClinChek® Controls are used for internal quality control in clinical chemical laboratories. The Hemoglobin F/A2/S-Control is a ready-to-use hemolysate of human blood with mean values in the elevated range.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After opening of the vials, the analytes are stable for 4 weeks, at 2 - 8 °C.

Notes:

Analysis of the control does not require hemolysis, that means you can inject the material directly into your analysing system. According to quality assurance all ClinChek® controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees exactly the same composition and stability for each vial of a lot.

Mean values:

The mean values and confidence intervals have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. The values for Hemoglobin F and A2 are traceable to the WHO-standards, NIBSC codes 85/616 and 89/666, respectively.

Pack size:

ClinChek® Hemolysate Control Level I
5 x 3 ml, **order no.: 13183**

Precautions:

The human whole blood which was used for manufacturing these control was tested for the absence of Hepatitis B-surface antigen and was HIV-negative. Nevertheless, the control should be considered as potentially infectious and treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH