ClinChek[®] - Control
Urine Control lyophilised / Kontrollurin lyophilisiert
FOR PORPHYRINS / FÜR PORPHYRINE

Intended use:

ClinChek[®] Urine Controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human urine and are available with mean values in the normal as well as in the pathological range. After reconstitution the controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 5.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

This product will be stable until the expiration date when stored unopened below -18 °C in the dark.

After reconstitution the stability of the analytes is:

- 1 day at 15 - 30 °C in the dark
- 3 days at 2 - 8 °C in the dark
- 90 days below -18 °C in the dark (avoid repeated freezing and thawing)

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek[®] Controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.68 %.

Mean values:

The mean values and confidence intervals have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinChek[®] Urine Control
Level I
10 x 5 ml, **order no.: 17080**
Level II
10 x 5 ml, **order no.: 17081**
Level I, II
2 x 5 x 5 ml, **order no.: 17082**

Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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- (1) RECIPE-HPLC
(2) RECIPE-Ion exchange chromatography / Photometry
RECIPE-Ionenaustausch-Chromatographie / Photometrie

RECIPE CHEMICALS + INSTRUMENTS GmbH
Dessauerstr. 3 D-80992 Munich / Germany
Tel.: +49 / 89 / 54 70 81 - 0 Fax: +49 / 89 / 54 70 81 - 11
Internet: www.recipe.de e-mail: info@recipe.de