## ClinChek® - Control Plasma Control lyophilised / Kontrollplasma lyophilisiert

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#### Intended use:

ClinChek® Plasma Controls are used for internal quality assurance in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. These lyophilised controls are based on human plasma and are available in two different ranges of concentration. After reconstitution the controls have to be prepared like patient samples in one series of analyses.

#### Reconstitution:

Add exactly 3.0 ml of analytical grade or double distilled water to the vial and close it well. Then dissolve the lyophilised material completely. Prior to use mix the controls well again.

#### Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C. After reconstitution the stability of the analytes is:

- 8 hours when stored at 15 - 30 °C
- 4 days when stored at 2 - 8 °C
- 30 days when stored below -18 °C (avoid repeated freezing and thawing)

#### Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing. RECIFE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %.

#### Mean values:

The mean values and confidence intervals have been established in independent reference laboratories with supervision of RECIFE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. For the determination of the trace elements, atomic absorption spectrometry (AAS) was used as well as ICP-MS. Even if you use different techniques of analysis, your results should be within the control ranges indicated on the data sheet.

#### Pack size:

ClinChek® Plasma Control  
Level I  
10 x 3 ml, **order no.: 8883**  
Level II  
10 x 3 ml, **order no.: 8884**  
Level I, II  
2 x 5 x 3 ml, **order no.: 8885**

#### Precautions:

The human plasma which was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen. Nevertheless, the plasma controls should be considered as potentially infectious and treated with appropriate care.

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In Germany Lithium and Magnesium are subject to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK). For assessment of the internal quality control, please also refer to table B 1a of the guideline. /

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RECIPE CHEMICALS + INSTRUMENTS GmbH