ClinChek® - Control Urine Control lyophilised / Kontrollurin lyophilisiert

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Intended use:

ClinChek® Urine Controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human urine and are available in three different ranges of concentration. After reconstitution the controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 1.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C.

After reconstitution the stability of the analytes is:

- 1 day at 15 - 30 °C
- 7 days at 2 - 8 °C
- 90 days below -18 °C
(avoid repeated freezing and thawing)

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® Controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %.

The creatinine concentration of the controls is: 0.22 g/l (level I, II and III).

Mean values:

The mean values and confidence intervals have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods.

Pack size:

ClinChek® Urine Control
Level I
10 x 1 ml, **order no.: MS8080**
Level II
10 x 1 ml, **order no.: MS8081**
Level III
10 x 1 ml, **order no.: MS8082**
Level I, II, III
3 x 3 x 1 ml, **order no.: MS8083**

Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH